

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT		MX 447 L		CONTROL		L		(Exp.)			2024-08-05			
		Rev 2		ABX Lyse										
PARAMETRES PARAMETERS	UNITES UNITS	MICROS			MICROS			PENTRA			TOLERANCES TOLERANCE			
		ES60 Care ST												
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0										± 0.4		
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	2.28										± 0.17		
HB HGB	g/dl	6.1										± 0.5		
	g/l	61										± 5		
	mmol/l	3.79										± 0.31		
HT HCT	%	16.4										± 2.0		
	l/l	0.164										± 0.020		
VGM MCV	µm ³ ; fl	72										± 4		
TGMH MCH	pg	26.8										± 2.0		
	fmol	1.66										± 0.12		
CCMH MCHC	g/dl	37.2										± 3.0		
	g/l	372										± 30		
	mmol/l	23.08										± 1.86		
IDR RDW	%	15.0										± 3.0		
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	79										± 22		
VPM MPV	µm ³ ; fl	8.4										± 2.0		
LYMPHO	#	1.22										± 0.40		
	%	61.0										± 8.0		
MONO	#	0.21										± 0.19		
	%	10.5										± 9.5		
GRANULO	#	0.57										± 0.40		
	%	28.5										± 10.0		

LOT		MX 447 N		CONTROL		N		(Exp.)			2024-08-05			
		Rev 2		ABX Lyse										
PARAMETRES PARAMETERS	UNITES UNITS	MICROS			MICROS			PENTRA			TOLERANCES TOLERANCE			
		ES60 Care ST												
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.5										± 0.8		
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	4.56										± 0.26		
HB HGB	g/dl	13.4										± 0.7		
	g/l	134										± 7		
	mmol/l	8.32										± 0.43		
HT HCT	%	36.9										± 2.7		
	l/l	0.369										± 0.027		
VGM MCV	µm ³ ; fl	81										± 4		
TGMH MCH	pg	29.4										± 2.0		
	fmol	1.82										± 0.12		
CCMH MCHC	g/dl	36.3										± 3.0		
	g/l	363										± 30		
	mmol/l	22.53										± 1.86		
IDR RDW	%	15.5										± 3.0		
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	264										± 40		
VPM MPV	µm ³ ; fl	7.5										± 2.0		
LYMPHO	#	2.63										± 0.40		
	%	35.0										± 5.0		
MONO	#	0.64										± 0.58		
	%	8.5										± 7.7		
GRANULO	#	4.24										± 0.75		
	%	56.5										± 10.0		

LOT		MX 447 H		CONTROL		H		(Exp.)			2024-08-05			
		Rev 2		ABX Lyse										
PARAMETRES PARAMETERS	UNITES UNITS	MICROS			MICROS			PENTRA			TOLERANCES TOLERANCE			
		ES60 Care ST												
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	20.4										± 1.7		
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	5.78										± 0.32		
HB HGB	g/dl	18.1										± 0.9		
	g/l	181										± 9		
	mmol/l	11.24										± 0.56		
HT HCT	%	51.4										± 3.4		
	l/l	0.514										± 0.034		
VGM MCV	µm ³ ; fl	89										± 4		
TGMH MCH	pg	31.3										± 2.0		
	fmol	1.94										± 0.12		
CCMH MCHC	g/dl	35.2										± 3.0		
	g/l	352										± 30		
	mmol/l	21.85										± 1.86		
IDR RDW	%	15.0										± 3.0		
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	509										± 65		
VPM MPV	µm ³ ; fl	7.4										± 2.0		
LYMPHO	#	3.06										± 1.00		
	%	15.0										± 5.0		
MONO	#	1.96										± 1.76		
	%	9.6										± 8.6		
GRANULO	#	15.38										± 2.04		
	%	75.4										± 10.0		