

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT		MX 447 L		CONTROL		L		(Exp.)		2024-08-05		
		Rev 2								(YYYY-MM-DD)		
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE		
		MICROS										
		EMI CRP										
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l									2.0		± 0.4
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l									2.29		± 0.15
HB HGB	g/dl									6.1		± 0.4
	g/l									61		± 4
	mmol/l									3.79		± 0.25
HT HCT	%									16.3		± 2.0
	l/l									0.163		± 0.020
VGM MCV	µm ³ ; fl									71		± 4
TGMH MCH	pg									26.6		± 2.0
	fmol									1.65		± 0.12
CCMH MCHC	g/dl									37.5		± 3.0
	g/l									375		± 30
	mmol/l									23.30		± 1.86
IDR RDW	%									14.5		± 3.0
PLAQ_PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l									74		± 20
VPM MPV	µm ³ ; fl									7.4		± 2.0
LYMPHO	#									1.23		± 0.40
	%									61.5		± 8.0
MONO	#									0.19		± 0.17
	%									9.5		± 8.6
GRANULO	#									0.58		± 0.40
	%									29.0		± 10.0

LOT		MX 447 N		CONTROL		N		(Exp.)		2024-08-05		
		Rev 2								(YYYY-MM-DD)		
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE		
		MICROS										
		EMI CRP										
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l									7.7		± 0.8
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l									4.50		± 0.18
HB HGB	g/dl									13.4		± 0.5
	g/l									134		± 5
	mmol/l									8.32		± 0.31
HT HCT	%									36.5		± 2.5
	l/l									0.365		± 0.025
VGM MCV	µm ³ ; fl									81		± 4
TGMH MCH	pg									29.8		± 2.0
	fmol									1.85		± 0.12
CCMH MCHC	g/dl									36.8		± 3.0
	g/l									368		± 30
	mmol/l									22.83		± 1.86
IDR RDW	%									14.0		± 3.0
PLAQ_PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l									257		± 40
VPM MPV	µm ³ ; fl									7.5		± 2.0
LYMPHO	#									2.70		± 0.40
	%									35.0		± 5.0
MONO	#									0.58		± 0.52
	%									7.5		± 6.8
GRANULO	#									4.43		± 0.77
	%									57.5		± 10.0

Ref:TEMP-0830 Rev.43 BACK/VEISO 1300104144

LOT		MX 447 H		CONTROL		H		(Exp.)		2024-08-05		
		Rev 2								(YYYY-MM-DD)		
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE		
		MICROS										
		EMI CRP										
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l									20.5		± 1.6
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l									5.67		± 0.20
HB HGB	g/dl									18.0		± 0.6
	g/l									180		± 6
	mmol/l									11.18		± 0.37
HT HCT	%									49.9		± 3.0
	l/l									0.499		± 0.030
VGM MCV	µm ³ ; fl									88		± 4
TGMH MCH	pg									31.7		± 2.0
	fmol									1.97		± 0.12
CCMH MCHC	g/dl									36.1		± 3.0
	g/l									361		± 30
	mmol/l									22.40		± 1.86
IDR RDW	%									13.5		± 3.0
PLAQ_PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l									481		± 45
VPM MPV	µm ³ ; fl									7.4		± 2.0
LYMPHO	#									3.38		± 1.00
	%									16.5		± 5.0
MONO	#									2.62		± 2.36
	%									12.8		± 11.5
GRANULO	#									14.49		± 2.05
	%									70.7		± 10.0