

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT		MX 447 L		CONTROL		L		(Exp.)		2024-08-05	
		Rev 2								(YYYY-MM-DD)	
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		ES60	MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	
							60 BP - 60 18P				
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0				2.0					± 0.4
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	2.23				2.26					± 0.17
HB HGB	g/dl	5.8				5.7					± 0.5
	g/l	58				57					± 5
	mmol/l	3.60				3.52					± 0.31
HT HCT	%	16.1				16.5					± 2.0
	l/l	0.161				0.165					± 0.020
VGM MCV	µm ³ ; fl	72				73					± 4
TGMH MCH	pg	26.0				25.1					± 2.0
	fmol	1.62				1.56					± 0.12
CCMH MCHC	g/dl	36.1				34.4					± 3.0
	g/l	361				344					± 30
	mmol/l	22.43				21.35					± 1.86
IDR RDW	%	15.0				15.0					± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	72				79					± 22
VPM MPV	µm ³ ; fl	8.4				7.4					± 2.0
LYMPHO	#	1.24				1.28					± 0.40
	%	62.0				64.0					± 8.0
MONO	#	0.21				0.21					± 0.19
	%	10.5				10.5					± 9.5
GRANULO	#	0.55				0.51					± 0.40
	%	27.5				25.5					± 10.0

LOT		MX 447 N		CONTROL		N		(Exp.)		2024-08-05	
		Rev 2								(YYYY-MM-DD)	
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		ES60	MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	
							60 BP - 60 18P				
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.5				7.3					± 0.8
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	4.51				4.53					± 0.26
HB HGB	g/dl	13.4				13.2					± 0.7
	g/l	134				132					± 7
	mmol/l	8.32				8.22					± 0.43
HT HCT	%	37.0				37.1					± 2.7
	l/l	0.370				0.371					± 0.027
VGM MCV	µm ³ ; fl	82				82					± 4
TGMH MCH	pg	29.7				29.2					± 2.0
	fmol	1.85				1.81					± 0.12
CCMH MCHC	g/dl	36.2				35.6					± 3.0
	g/l	362				356					± 30
	mmol/l	22.50				22.12					± 1.86
IDR RDW	%	15.0				15.0					± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	244				259					± 40
VPM MPV	µm ³ ; fl	7.8				6.7					± 2.0
LYMPHO	#	2.59				2.63					± 0.40
	%	34.5				36.0					± 5.0
MONO	#	0.64				0.62					± 0.56
	%	8.5				8.5					± 7.7
GRANULO	#	4.28				4.05					± 0.75
	%	57.0				55.5					± 10.0

LOT		MX 447 H		CONTROL		H		(Exp.)		2024-08-05	
		Rev 2								(YYYY-MM-DD)	
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		ES60	MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	
							60 BP - 60 18P				
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	20.3				19.7					± 1.7
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	5.76				5.70					± 0.32
HB HGB	g/dl	18.6				18.1					± 0.9
	g/l	186				181					± 9
	mmol/l	11.55				11.25					± 0.56
HT HCT	%	51.3				51.3					± 3.4
	l/l	0.513				0.513					± 0.034
VGM MCV	µm ³ ; fl	89				90					± 4
TGMH MCH	pg	32.3				31.8					± 2.0
	fmol	2.01				1.97					± 0.12
CCMH MCHC	g/dl	36.3				35.3					± 3.0
	g/l	363				353					± 30
	mmol/l	22.53				21.93					± 1.86
IDR RDW	%	14.5				14.5					± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	478				494					± 65
VPM MPV	µm ³ ; fl	7.6				6.6					± 2.0
LYMPHO	#	3.25				3.05					± 1.00
	%	16.0				15.5					± 5.0
MONO	#	1.95				1.90					± 1.71
	%	9.6				9.6					± 8.6
GRANULO	#	15.10				14.75					± 1.97
	%	74.4				74.9					± 10.0